

**SPINE —
CORRECTION
SOLUTION**




spinamic

CE FDA Pmda

#1302, 13F Flag One
Nambusunhwan-ro 2621
Gangnam-gu, Seoul
Phone. 82-2-6447-9124
Fax. 82-2-6442-9124
E-mail. contact@vntc.me
www.vntc.me


PATIENT CASE 1



Non-brace
Cobb angle 28°

spinamic®
Cobb angle 14°

PATIENT CASE 2



Non-brace
Cobb angle 22°

spinamic®
Cobb angle 11°

HARD BRACE	Maximized correction
	<ul style="list-style-type: none"> ✗ Pain control ✗ Wearability ✗ Mobility ✗ Adjustability ✗ Washability
	Constricts bone growth & breathing, causes rib cage damage, muscle weakness, and mental depression

PATIENT'S GOAL IMPROVED

SPINAMIC	Accurate and smart correction program based on individual patient data analysis
	<ul style="list-style-type: none"> ✓ Pain reduction ✓ Wearability ✓ Mobility ✓ Adjustability ✓ Washability
	Minimized side effects Patients' quality of life Improved


FEATURE #1



THREE POINT PRESSURE

Spinamic brace generates desirable mechanical force based on the 3 points lateral compression principle to achieve the whole body alignment and balance.


FEATURE #2



COMPRESSION DIALS

On both sides of **Spinamic** brace, there are compression-strength-adjustable dials to provide appropriate amount of therapy for each patient.

FEATURE #3



SPINAMIC APPLICATION

We apply the most advanced solution to improve patient compliance and get the clinical output we want. **Spinamic** app provides push notification function to help the patients keep track of the device wearing time.

**BEAUTIFUL,
COMFORTABLE,
AND ALSO
FUNCTIONAL
TEXTILES
FOR IMPROVED
MOVEMENT.**

